

# Farnham Lab Protocol for CRISPR/Cas 9-mediated enhancer deletion using a puromycin selection vector

Created by Yu (Phoebe) Guo 20151024

## Part A. Design gRNA oligos

1. Design gRNAs  
Go to Optimized CRISPR design website: <http://crispr.mit.edu/>. Input 23-500nt genomic sequence, choose the correct target genome and submit the job. Choose gRNA that has a score greater than 85.
2. gRNA specificity check
  - a) Go to CRISPR RGEN tool Cas-OFFinder website: <http://www.rgenome.net/cas-offinder/>. Select PAM type (5'-NGG-3') and target genome. Input the 20nt gRNA w/o PAM as query sequence. Allow mismatch number equal or less than 2. Submit the job and check result. Avoid using a gRNA that has 1 or 2 mismatch sequences.
  - b) Back to Cas-OFFinder again, input the 15nt gRNA that close to PAM (w/o PAM) as query sequence. Allow mismatch number equal or less than 2. Submit the job and check result. Avoid using gRNA has other 0 mismatch sequences since the sequence near PAM is critical for CRISPR specificity.
3. Order oligos  
For pSpCas9(BB)-2A-Puro (PX459) V2.0 (*Addgene plasmid #62988*) backbone vector, order reverse complement oligos that contains 5' and 3' overhangs as shown below:  
e.g.  
20nt gRNA sequence: TGTTCTTCTTATCGTGCCAT  
Forward oligo: CACCGTGTTCCTTCTTATCGTGCCAT  
Reverse oligo: AAACATGGCACGATAAGAAGAACAC

## Part B. Plasmid cloning

1. Annealing of oligo

### Option A

- a) In a PCR tube, mix:

100µM Forward oligo	15 µl
100µM Reverse oligo	15 µl
Annealing buffer	45 µl
<hr/>	
Total Reaction Volume	75 µl

\*\*Annealing buffer (10mM Tris-HCl pH8.0, 1mM EDTA, 100mM NaCl).

- b) Incubate in thermal cycler 92°C 2min, 72°C 2min, 55°C 2min, 37°C 2min, hold at 10°C.

OR

**Option B**

a) In a PCR tube, mix:

100µM Forward oligo	1 µl
100µM Reverse oligo	1 µl
10X T4 Ligation Buffer (NEB)	1 µl
ddH <sub>2</sub> O	6.5 µl
T4 PNK	0.5 µl
<hr/>	
Total Reaction Volume	10 µl

b) Incubate in thermal cycler 37°C 30min, 95°C 5min, ramp down to 25°C at 5°C/min, hold at 10°C.

2. Digest pSpCas9(BB)-2A-Puro (PX459) V2.0 (*Addgene plasmid #62988*) with BbsI endonuclease (*NEB Cat# R0539S/R0539L*) or FastDigest BbsI (*Thermo Cat.#FD1014*), reaction system:

Restriction Enzyme	NEB BbsI 1ul ( or 10 units)	FastDigest BbsI 1ul
PX459 Vector	1 µg	2ul (up to 1 µg)
10X Buffer	5 µl	2 µl
ddH <sub>2</sub> O	x µl	x µl
<hr/>		
Total Reaction Volume	50 µl	20 µl

Incubate at 37°C for 1-2 hrs for NEB enzyme, 5-30 mins for FastDigest enzyme.

3. Gel purification of BbsI-digested PX459 vector using QIAquick Gel extraction kit (*QIAGEN, Cat# 28704/28706*).

4. Ligate annealed oligo duplex with digested PX459 vector  
 a. If using **Option A**, use T4 ligase (*NEB #M0202*) in this step.

Diluted oligo duplex from step 1 (1:20)	1 µl
Digested PX459 vector	50 ng
10X T4 ligase buffer (NEB)	2 µl
T4 ligase	1 µl
ddH <sub>2</sub> O	x µl
<hr/>	
Total Reaction Volume	20 µl

Incubate at 16°C overnight

If using **Option B**, use Quick Ligase (*NEB #M2200*) in this step:

Diluted oligo duplex from step 1 (1:250)	1 µl
Digested PX459 vector	50 ng
2X Quick Ligation buffer (NEB)	5 µl
Quick ligase	1 µl
ddH <sub>2</sub> O	x µl
<hr/>	
Total Reaction Volume	11 µl

Incubate at 37°C for 10 min

5. OPTIONAL: Treat ligation reaction with PlasmidSafe exonuclease to digest any residual linearize DNA:

Ligation product from previous step	11 $\mu$ l
10X PlasmidSafe Buffer	1.5 $\mu$ l
10Mm ATP	1.5 $\mu$ l
PlasmidSafe exonuclease	1 $\mu$ l
Total Reaction Volume	15 $\mu$ l

Incubate at 37 °C for 30 min, followed by 70 °C for 30 min. After PlasmidSafe treatment, the reaction can be stored at -20 °C for at least 1 week.

6. Transformation of PX459 ligation product by using XL-10 gold Ultracompetent cells (*Agilent Technologies, Cat# 200315*)
  - a) Pre-chill several 14ml BD Falcon polypropylene round-bottom tubes on ice. (Remember to include one tube for negative control and one tube for the pUC18 positive control.) Preheat SOC medium to 42°C.
  - b) Thaw the competent cells on ice. When thawed, gently mix and aliquot 100  $\mu$ l of cells into each of the pre-chilled tubes.
  - c) Add 4  $\mu$ l of the  $\beta$ -ME mix provided with this kit to each aliquot of cells.
  - d) Swirl the tubes gently. Incubate the cells on ice for 10 minutes, swirling gently every 2 minutes.
  - e) Add 0.1-50 ng of the experimental DNA (or 2  $\mu$ l of a ligation mixture) to one aliquot of cells. Dilute the pUC18 control DNA (provided in the kit) 1:10 with sterile dH<sub>2</sub>O, then add 1  $\mu$ l of the diluted pUC18 DNA to the other aliquot of cells as positive control. Add 2  $\mu$ l of sterile dH<sub>2</sub>O to negative control
  - f) Swirl the tubes gently, then incubate the tubes on ice for 30 minutes.
  - g) Heat-pulse the tubes in a 42°C water bath for 30 seconds. The duration of the heat pulse is critical.
  - h) Incubate the tubes on ice for 2 minutes.
  - i) Add 0.9 ml of preheated (42°C) SOC medium and incubate the tubes at 37°C for 1 hour with shaking at 225-250 rpm.
  - j) Plate  $\leq$ 200  $\mu$ l of the transformation mixture on LB agar plates containing the appropriate antibiotic (in this case, ampicillin). For the pUC18 control transformation, plate 5  $\mu$ l of the transformation on LB-ampicillin agar plates.
  - k) Incubate the plates at 37°C overnight.
  - l) For the pUC18 control, expect 250 colonies ( $\geq 5 \times 10^9$  cfu/  $\mu$ g pUC18 DNA). For the experimental DNA, the number of colonies will vary according to the size and form of the transforming DNA, with larger and non-supercoiled DNA producing fewer colonies.
7. Isolation of recombinant plasmid by using QIAprep Spin Miniprep Kit(250) (*Qiagen, Cat# 27106*).
8. Send plasmid for Sanger sequencing to confirm successful cloning.

### Part C. Isolation of CRISPR Knock out clonal cell line

1. Transfection of recombinant plasmid to cells using Lipofectamine LTX reagent (*Invitrogene, Cat# 15388-100*) or other suitable transfection reagents.
  - a) Plate cells ( $3.0 \times 10^5$ /well) in 6well plate 1-2 days to reach 60-70% confluence.

- \*\*Make sure you have one well for no gRNA negative control PX459 plasmid.
- \*\*GFP plasmid can be used to test transfection efficiency.
- b) One the day of transfection, in a tube dilute Lipofectamine LTX reagent 5-12.5ul/well in Opti-MEM reduce serum medium (*ThermoFisher, Cat.# 31985062*) 125ul/well.
- \*\*5ul/well works for C4-2B and LNCaP cells
- c) In a tube, dilute recombinant plasmid(s) or negative control PX459 plasmid 2500ng/well in Opti-MEM medium 125 ul/well, add Plus Reagent 2.5 ul/well.
- \*\*If using multiple gRNAs, add plasmids at a 1:1 ratio to a total DNA amount of 2500ng/well
- d) Add diluted DNA to diluted Lipofectamine LTX reagent at 1:1 ratio (125ul+125ul/well).
- e) Incubate 5min at RT to make DNA-lipid complex.
- f) Take out old medium from the wells
- g) Add DNA-lipid complex 250ul/well to cells with 2ml/well fresh complete growth medium
- h) Incubate at 37°C overnight
- 2. Gently replace medium with 1ml of fresh complete growth medium + puromycin
  - \*\*For C4-2B cells, use 400-900ng/ml of puromycin
- 3. Incubate at 37°C for 48-72hrs, until most of cells in un-transfected well are dead
- 4. Disassociate the transfected cells to single cell suspension
- 5. Count cells and serially dilute the cells to 0.3-0.5 cell/ 100ul (that's 3-5 cells/ml), plate the single cell suspension into several 96well plates (I usually do 2-3 plates) 100ul/well, return to incubator.
  - \*\*Conditional medium (40% old growth medium + 60% complete medium) might be needed for some cell lines at the last dilution.
  - e.g. RWPE-1, RWPE-2

## OR

Sort live cells into 96 well plates contain 100ul/well growth medium at 1cell/well, return to incubator.

6. Freeze down the leftover transfected cell population for future potential use.
7. OPTIONAL: Test deletion efficiency on the population using leftover cells.
  - a) Dissociate transfected cell population and spin them down at 200g for 5 min at room temperature.
  - b) Aspirate the medium completely
  - c) Extract DNA by adding QuickExtract DNA extraction solution (*Epicentre, cat.#QE09050*), 50 µl or 10 µl of the solution for each well of a 24-well or 96-well plate, respectively, is typically used.
  - d) Vortex mix for 15 seconds.
  - e) Transfer the tube to 65°C and incubate for 6 minutes.
  - f) Vortex mix for 15 seconds.
  - g) Transfer the tube to 98°C and incubate for 2 minutes.
  - h) Store the DNA at -20°C, or at -70°C for long term storage.
  - i) PCR using primers flank the deletion region, use 5 µl or less of the extracted DNA for each PCR amplification (1 µl is usually used).

8. Allow several weeks (2-3 weeks usually) for each single cell to grow into a visible colony, inspect cells every week and change medium when necessary.
9. Disassociate the colonies and re-plate half of them into 96well plates to allow cells to spread out and avoid clumping.
10. Harvest other half of cells to test for deletion
  - \*\*Make sure to label the tubes correspondingly to the wells
  - a) Spin them down at 200g for 5 min at room temperature.
  - b) Aspirate the medium completely.
  - c) Extract DNA by adding 10  $\mu$ l QuickExtract DNA extraction.
  - d) Vortex mix for 15 seconds.
  - e) Transfer the tube to 65°C and incubate for 6 minutes.
  - f) Vortex mix for 15 seconds.
  - g) Transfer the tube to 98°C and incubate for 2 minutes.
  - h) Store the DNA at -20°C, or at -70°C for long term storage.
  - i) PCR using primers flank the deletion region as well as primers inside of deletion region, use 5  $\mu$ l or less of the extracted DNA for each PCR amplification (1  $\mu$ l is usually used).
  - j) For colonies that show no inside band and a clean deletion band, congratulations, you just got yourself a complete deletion clone!
11. Keep growing and expanding the deletion clone for storage and future functional assays.

#### **REFERENCE**

Ran, F. Ann, et al. "Genome engineering using the CRISPR-Cas9 system." *Nature protocols* 8.11 (2013): 2281-2308.